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
Levenhuk 1ST
Levenhuk 2ST
Levenhuk 3ST
Levenhuk 4ST

- EN User Manual
- BG Ръководство за потребителя
- CZ Návod k použití
- DE Bedienungsanleitung
- ES Guía del usuario
- HU Használati útmutató
- IT Guida all'utilizzo
- PL Instrukcja obsługi
- PT Manual do usuário
- RU Инструкция по эксплуатации


*Наслади се отблизо
Radost zaostřit
Zoom ran und hab Fun!
Amplie y disfrute
Kellemes nagyítást!
Ingrandisci il divertimento
Radość przybliżania
Dê um zoom na sua emoção
Приближает с удовольствием*


levenhuk
Zoom&Joy


1a. Levenhuk 1ST

1b. Levenhuk 2ST

1c. Levenhuk 3ST

1d. Levenhuk 4ST

EN

1. Pole
2. Locking knob
3. Focus knob
4. Eyepieces
5. Diopter adjustment ring
6. Binocular head
7. Objective(s)
8. Stage clips
9. Stage
10. Upper illumination (reflected light)
11. Battery compartment
12. Illumination switch
13. Glass stage inset
14. Lower illumination (transmitted light, in the base)
15. Brightness adjustment wheel
16. Power connector
17. Power switch

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Levenhuk ST Series Stereo Microscopes

Please read this manual carefully before using your microscope for the first time.

EN

General Information

Levenhuk 1ST, Levenhuk 2ST, Levenhuk 3ST and Levenhuk 4ST microscopes were designed to be used during assembly of intricate equipment, including electronics; for school experiments and scientific research. All models may be found in schools and research institutes, at home and in manufacturing plants.

Using the microscope

Illumination (Levenhuk 1ST, Levenhuk 3ST, Levenhuk 4ST)

Plug your microscope into a power source or insert batteries, turn it on and choose the required illumination source. If your microscope is fitted with brightness control, you can also choose the brightness level of illumination.

Stage inset glass

- Frosted glass stage inset may be used for observations of transparent specimens or thin slices of objects. Place it on the stage and secure it in place with the locking screw. Lower illumination source should be used during observations with this glass stage inset.
- The kit of your new microscope may also include a black and white glass stage inset. Normally, you would want to use the white side; however, if a specimen is very bright, or colored white, you should use the black side for additional contrast. Upper illumination source should be used during observations with this glass stage inset. Before installing the black and white glass make sure to remove the frosted glass from the stage.

Preparations

Place a specimen in the center of a blank slide, put the slide on the stage and lock it in place with stage clips.

Eyecups

For additional comfort during observations you can use rubber eyecups on both eyepieces.

Focusing, diopter adjustment, interpupillary distance

Place a slide on the stage and lock it in place with stage clips. Loosen the locking knob on the microscope pole and slide the head up or down into a desired position and retighten the locking knob. Look through the right eyepiece and turn the objective lens handle until you see the specimen in your field of view. Afterwards, turn the focus knob until the image is sharp and clear. Look through the left eyepiece and rotate the diopter adjustment ring until the image is sharp and clear. Do not adjust focus during diopter adjustment! Upon completing this process, look through both eyepieces, hold them with your hands and bring them closer or move them apart until a double image turns into a single image.

Changing objective lenses

Levenhuk 3ST and Levenhuk 4ST has a revolving nosepiece with two objective lenses. To switch between them, you have to rotate the nosepiece.

Changing lamps

CAUTION! Before changing a lamp or a fuse make sure the microscope is unplugged from a power source. Also, wait until the lamp cools down.

Changing the upper illumination lamp

Loosen the locking screw on the upper illumination housing and remove the light socket from within. Replace the old lamp with a new one and put the light socket back. Retighten the locking screw.

Changing the lower illumination lamp

Loosen the locking screw on the stage and remove the glass. Remove the old lamp from the light socket and replace it with a new one. Put the glass back in its place and retighten the locking screw.

Changing the fuse replacement

The fuse box is located on the bottom of the microscope base. Remove the cover and replace the old fuse with a new fuse.

Specifications

	Levenhuk 1ST	Levenhuk 2ST	Levenhuk 3ST	Levenhuk 4ST
Head	binocular, vertical		binocular, inclined at 45°	
Optics material	optical glass			
Magnification	20x	40x	20–40x	
Eyepiece tube diameter	23.2mm	30.5mm		30mm
Eyepieces	WF10x (2pcs)			
Objectives	2x	4x	2x, 4x	
Working distance	65mm	60mm	57mm	130mm
Interpupillary distance	55–75mm	60mm	60mm	55–75mm
Glass stage inset	50mm, with stage clips, white stage plate	50mm, with stage clips, double-sided (black and white, plastic) stage plate	90mm, matte glass stage with stage clips, additional double-sided (black and white, plastic) stage plate	95mm, matte glass stage with stage clips, additional double-sided (black and white, plastic) stage plate
Eyepiece diopter adjustment	–	±5 D		
Focus	coarse, 35mm	coarse, 40mm		coarse, 35mm
Power supply	2pcs AA batteries (not included in the kit)	–	220V/50Hz or 110V/60Hz (optional)	
Illumination	upper (0.06W LED)	natural illumination	upper and lower halogen, 10W	upper (1W LED) lower (2W LED)
Operating temperature range	–5...40 °C (23...104 °F)			5...35 °C (41...95 °F)
Additional accessories	dust cover	rubber eyecups (2pcs)	rubber eyecups (2pcs), backup halogen lamp 12V/10W	rubber eyecups (2pcs), dust cover

The manufacturer reserves the right to make changes to the product range and specifications without prior notice.

Caution: Please refer to the specifications table for the correct mains voltage and never attempt to plug a 110V device into 220V outlet and vice versa without using a converter. Remember that mains voltage in the U.S. and Canada is 110V and 220–240V in most European countries.

Care and maintenance

- **Never, under any circumstances, look directly at the Sun, another bright source of light or at a laser through this device, as this may cause PERMANENT RETINAL DAMAGE and may lead to BLINDNESS.**
- Take necessary precautions when using the device with children or others who have not read or who do not fully understand these instructions.
- After unpacking your microscope and before using it for the first time check for integrity and durability of every component and connection.
- Do not try to disassemble the device on your own for any reason. For repairs and cleaning of any kind, please contact your local specialized service center.
- Protect the device from sudden impact and excessive mechanical force. Do not apply excessive pressure when adjusting focus. Do not overtighten the locking screws.
- Do not touch the optical surfaces with your fingers. To clean the device exterior, use only special cleaning wipes and special optics cleaning tools from Levenhuk. Do not use any corrosive or acetone-based fluids to clean the optics.
- Abrasive particles, such as sand, should not be wiped off lenses, but instead blown off or brushed away with a soft brush.
- Do not use the device for lengthy periods of time, or leave it unattended in direct sunlight. Keep the device away from water and high humidity.
- Be careful during your observations, always replace the dust cover after you are finished with observations to protect the device from dust and stains.
- If you are not using your microscope for extended periods of time, store the objective lenses and eyepieces separately from the microscope.
- Store the device in a dry, cool place away from hazardous acids and other chemicals, away from heaters, open fire and other sources of high temperatures.
- When using the microscope, try not to use it near flammable materials or substances (benzene, paper, cardboard, plastic, etc.), as the base may heat up during use, and might become a fire hazard.
- Always unplug the microscope from a power source before opening the base or changing the illumination lamp. Regardless of the lamp type (halogen or incandescent), give it some time to cool down before trying to change it, and always change it to a lamp of the same type.
- Always use the power supply with the proper voltage, i.e. indicated in the specifications of your new microscope. Plugging the instrument into a different power outlet may damage the electric circuitry of the microscope, burn out the lamp, or even cause a short circuit.
- **Seek medical advice immediately if a small part or a battery is swallowed.**

Battery safety instructions

- Always purchase the correct size and grade of battery most suitable for the intended use.
- Always replace the whole set of batteries at one time; taking care not to mix old and new ones, or batteries of different types.
- Clean the battery contacts and also those of the device prior to battery installation.
- Make sure the batteries are installed correctly with regard to polarity (+ and –).
- Remove batteries from equipment that is not to be used for an extended period of time.
- Remove used batteries promptly.
- Never short-circuit batteries as this may lead to high temperatures, leakage, or explosion.
- Never heat batteries in order to revive them.
- Remember to switch off devices after use.
- Keep batteries out of the reach of children, to avoid risk of ingestion, suffocation, or poisoning.
- Do not disassemble batteries.
- Utilize used batteries as prescribed by your country laws.

Levenhuk International Lifetime Warranty

All Levenhuk telescopes, microscopes, binoculars and other optical products, except for accessories, carry a **lifetime warranty** against defects in materials and workmanship. Lifetime warranty is a guarantee on the lifetime of the product on the market. All Levenhuk accessories are warranted to be free of defects in materials and workmanship for **six months** from date of retail purchase. The warranty entitles you to free repair or replacement of the Levenhuk product in any country where a Levenhuk office is located if all warranty conditions are met.

For further details please visit our web site: www.levenhuk.com/warranty

If warranty problems arise, or if you need assistance in using your product, contact the local Levenhuk branch.

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Осветление (Levenhuk 1ST, Levenhuk 3ST, Levenhuk 4ST)

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